

**Protocol for Counting BMMs**  
**Cindy Rohde 1/7/04**

1. Aspirate media from well using a 200  $\mu$ L pipet tip.
2. Wash cells 3x with 1 mL ice-cold PBS.
3. Add appropriate volume (~300  $\mu$ L for a 24-well plate) of ice-cold PBS/0.005% Zwittergent to each well.
4. Incubate plate at 4°C for 10-15 min.
5. Use a 1 mL pipetman to detach cells from plate.
6. Check that most cells are detached under the microscope.
7. Use a 1 mL pipetman to resuspend cells and remove a 100  $\mu$ L aliquot.
8. Mix 100  $\mu$ L aliquot with 20  $\mu$ L trypan blue and transfer to hemacytometer.
9. Count cells.